

VISCOELASTIC

Synergy visc

Synergy visc Sodium Hyaluronate is a viscoelastic solution containing cohesive sodium hyaluronate for injection during eye surgery

- High biocompatibility
- Optical transparency
- Optimal stability in the anterior chamber and in capsular bag
- Perfect adhesion and protection of endothelial cells from mechanical damages
- Medium viscosity
- No refrigeration required
- Supplied in 1.1 ml syringe in kit with visco cannula



	1.0%	1.2%	1.4%	1.6%	1.8%	2.0%
pH	6.8 – 7.6	6.8 – 7.6	6.8 – 7.6	6.8 – 7.6	6.8 – 7.6	6.8 – 7.6
Osmolality	300 – 350 mOsm/Kg	300 – 350 mOsm/Kg	300 – 350 mOsm/Kg	300 – 350 mOsm/Kg	300 – 350 mOsm/Kg	300 – 350 mOsm/Kg
Dynamic viscosity	25.000 – 50.000 mPa.s	50.000 – 100.000 mPa.s	100.000 – 150.000 mPa.s	150.000 – 200.000 mPa.s	200.000 – 300.000 mPa.s	> 400.000 mPa.s
Molecular weight (Sodium Hyaluronate raw material)	2.5 x10 ⁶ - 3x10 ⁶ Da	2.5 x10 ⁶ - 3x10 ⁶ Da	2.5 x10 ⁶ - 3x10 ⁶ Da	2.5 x10 ⁶ - 3x10 ⁶ Da	2.5 x10 ⁶ - 3x10 ⁶ Da	2.5 x10 ⁶ - 3x10 ⁶ Da

Scan the QR code and download product technical data sheet

